

synchrony

Easy HD

Progressive Lenses

Easy to fit.
Easy to adapt.
Easy to view.



Understanding wearer's needs and wants.



66% of prebyopes often experience problems when using progressive lenses...

"I don't like it when I need to move my head to see things clearly."

"Objects seem to be moving when I suddenly look in different directions."

"It's difficult for me to adapt to new progressives."

A vast majority of presbyopes often experience problems with progressive lenses.

- ▶ Patients cannot see clearly at all distances, and objects look distorted.
- ▶ Reading or working at the computer can cause eye strain and fatigue.
- ▶ Patients experience blurred vision in the periphery, causing unnecessary head movements in order to see clearly.

Source: Market research run in 2014 with +1,000 users.

Easy to fit. Easy to adapt. Easy to view. Freeform lenses for all daily activities.

The synchrony Easy HD family of freeform progressive lenses provides patients who are looking for a value-for-money solution with multiple options:

Easy View HD:

- ▶ Large distortion-free distance vision.
- ▶ Enhanced near zone for easy reading.
- ▶ Available in 3 fitting height options: 17 mm, 15 mm and 13 mm.

Easy Adapt HD:

Minimum fitting height of 18 mm.

- ▶ Ultra-soft periphery for easy adaptation.
- ▶ Effective balance of central and peripheral viewing zones.
- ▶ Binocular optimization of the design.

Availability

Material	Diameter	Colour	Rx Range	Cyl	Add Power
1.50 Hard Resin	75/95	Clear	-10.00 to +6.00D	-6.00D	0.75 to 3.50D
1.53 Trivex®	72/92	Clear	-9.00 to +5.00D	-4.00D	0.75 to 3.50D
1.59 Polycarbonate	72/92	Clear	-9.00 to +6.00D	-6.00D	0.75 to 3.50D
1.60 High Index	75/95*	Clear	-12.00 to +6.00D	-6.00D	0.75 to 3.50D
1.67 High Index	70/90	Clear	-12.00 to +8.00D	-6.00D	0.75 to 3.50D
1.74 High Index	Varies	Clear	-16.00 to +10.00D	-6.00D	0.75 to 3.50D
PhotoFusion Self-Tinting Lenses					
1.50 PhotoFusion	75/95	Grey & Brown	-10.00 to +6.00D	-6.00D	0.75 to 3.50D
1.59 PhotoFusion	72/92	Grey & Brown	-9.00 to +6.00D	-6.00D	0.75 to 3.50D
1.60 PhotoFusion	75/95*	Grey & Brown	-12.00 to +6.00D	-6.00D	0.75 to 3.50D
1.67 PhotoFusion	70/90	Grey & Brown	-12.00 to +8.00D	-6.00D	0.75 to 3.50D
Transitions					
1.50 Transitions	75/95	Grey, Brown & Vantage	-10.00 to +6.00D	-6.00D	0.75 to 3.50D
1.53 Transitions	72/92	Grey & Brown	-9.00 to +5.00D	-4.00D	0.75 to 3.50D
1.59 Transitions	72/92	Grey, Brown & XTRActive	-9.00 to +6.00D	-6.00D	0.75 to 3.50D
1.60 Transitions	75/95*	Grey & Brown	-12.00 to +6.00D	-6.00D	0.75 to 3.50D
1.67 Transitions	70/90	Grey, Brown & XTRActive	-12.00 to +8.00D	-4.00D	0.75 to 3.50D
Sun & Polarized					
1.50 Polarized	75/95	Grey & Brown	-8.00 to +6.00D	-6.00D	0.75 to 3.50D
1.53 Trivex NXT® Polarized	72/92	Grey & Brown	-5.00 to +5.00D	-4.00D	0.75 to 3.50D
1.53 Trivex NXT Photo Polarized	72/92	Grey & Brown	-5.00 to +5.00D	-4.00D	0.75 to 3.50D
1.53 Trivex NXT Tint	72/92	NXT® sun tints, mirrors	-9.00 to +5.00D	-4.00D	0.75 to 3.50D
1.53 Trivex NXT Photochromic	72/92	NXT® sun colours	-9.00 to +5.00D	-4.00D	0.75 to 3.50D
1.59 Polarized	72/92	Grey & Brown	-9.00 to +6.00D	-6.00D	0.75 to 3.50D
1.60 Polarized	75/95*	Grey, Brown & G-15	-11.00 to +6.00D	-6.00D	0.75 to 3.50D
1.67 Polarized	70/90	Grey & Brown	-11.00 to +6.00D	-6.00D	0.75 to 3.50D

* Please confirm diameter availability for Rx's over +4.00D / -5.00D with your lab



synchrony

synchrony lenses offer
more designs,
more practice solutions,
more satisfied patients.

For more information contact your Carl Zeiss Vision representative
or visit www.synchronylenses.ca

German Precision Optics Manufactured in Canada.

Carl Zeiss Vision Inc.
1-800-268-6489